

70% IPA Wipes

General purpose disposable 70% Isopropyl Alcohol (IPA) wipes are ideal for a variety of cleaning applications, including electronics, tools, workspaces, and other surfaces that require quick removal of contaminants. The 18 cm x 14 cm (7" x 5.5") perforated wipes are constructed of durable, lint-free, tear-resistant polyester. Suitable for use in Class 6-7 Cleanrooms.

Wipes are presaturated in a solution of 70% Isopropyl Alcohol (CAS: 63-67-0) and 30% Deionized Water. The screw-top storage canister includes a convenient flip lid with pull-through dispenser. Each canister contains 100 wipes. Canister is foil-sealed to keep wipes fresh before first use.

Physical and Chemical Properties

Physical state:	Liquid
Odor:	Alcohol-like
Fabric material:	Polyester-Viscose
Fabric specifications:	45 gsm cross-spun lace (non-woven)
Appearance:	Clear, colorless liquid
Evaporation rate:	Slow
Volatility:	100%
Flashpoint and method:	20°C, Open Cup
Solubility in water:	Soluble
Shelf life:	2 years



Screw-top canister includes a convenient flip lid with pull-through dispenser.



Item: ISOWP18

Pack: Canister with 100 each of 18 cm x 14cm (7" x 5.5") perforated 70% IPA wipes

Case: 24 canisters per case



Roll contains 100 perforated polyester wipes.



Canister is foil-sealed to keep wipes fresh.



www.unitedsci.com

United Scientific Supplies, Inc.

3055 N. Oak Grove Ave. P: +1 847-336-7556
Waukegan IL 60087 F: +1 847-336-7571
www.unitedsci.com E: info@unitedsci.com