

3-Ply Ear Loop Facemask

Safe


- Bacterial Filtration Efficiency $\geq 95\%$
- Three layers of protection
- 25 gsm PP outer layer (blue)
- 22 gsm non-woven melt-blown filter layer
- 20 gsm PP inner layer (white)
- Non-sterile
- Produced in an FDA registered factory

Comfortable

- 3-ply lightweight design
- Soft ear loops
- Adjustable nose cover

Easy

- Disposable
- Flat, pleated design
- Simple to put on and take off



Item: MASK1-3P

Dimensions: 17.5 cm x 9.5 cm

Packing: Box of 50, 2000 per case



www.unitedsci.com

Email us at info@unitedsci.com
for more information.

United Scientific Supplies, Inc.

3055 N. Oak Grove Ave. P: +1 847-336-7556

Waukegan IL 60087

www.unitedsci.com

F: +1 847-336-7571

E: info@unitedsci.com