

KN95 Respirator Mask

Safe


- Particle Filtration Efficiency $\geq 95\%$
- Bacterial Filtration Efficiency $\geq 99\%$
- Five layers of protection
- Two non-woven fabric layers
- Two meltblown fabric layers
- One hot air cotton layer
- Non-sterile
- Latex-free
- Produced in an FDA registered factory

Comfortable

- Curved shape for easy breathing
- Soft ear loops
- Adjustable nose bridge strip

Easy

- Disposable
- Lightweight design
- Simple to put on and take off



Five layers of protection



Item: MASK2-K95

Dimensions: 15 cm x 21 cm

Packing: Box of 10



www.unitedsci.com

Email us at info@unitedsci.com
for more information.

United Scientific Supplies, Inc.

3055 N. Oak Grove Ave. P: +1 847-336-7556
Waukegan IL 60087 F: +1 847-336-7571
www.unitedsci.com E: info@unitedsci.com